

## NUCLEAR EXTRACTION (from cells)

1mM DTT

1X PIC (Protease Inhibitor Cocktail)

10uM Lactacystin

Extraction Buffer

Lysis Buffer

PBS

- 1) Add DTT, PIC, and lactacystin to lysis buffer.
- 2) Aspirate media.
- 3) Rinse cell with ice cold PBS (1x)
- 4) Add 1mL of lysis buffer / 10cm plate.
- 5) Scrape cells.
- 6) Pipette cell lysates into 1.5mL eppendorf tube.
- 7) Spin at 14K rpm for 1min at 4 °C.
- 8) Collect pellet (nuclear extract) or Supernatant (cytoplasmic extract) as desired.
- 9) Resuspend cells in 100ul of ice cold lysis buffer - mix by pipetting up and down.
- 10) Incubate on ice for 10min.
- 11) Spin cells @14K rpm for 10min. in a microfuge.
- 12) Separate supernatant (cyto) from pellet (nucleus).
- 13) Add DTT and PIC, and lactacystin to extraction buffer.
- 14) Add 80ul of extraction buffer into the pellet.
- 15) Resuspend by mixing up and down with pipette.
- 16) Shake @ slow pace for 10min. so the liquid goes up and down using the shaker in Welch's deli.
- 17) Spin @ 14K rpm for 10min. using a microfuge.
- 18) Separate pellet (DNA) from supernatant (nuclear protein).
- 19) **Supernatant is now your nuclear lysate!!**
- 20) Flash freeze lysate and store at -70 °C.

### NUCLEAR EXTRACTION BUFFER RECIPES

Lysis Buffer (200ml total)

1M Hepes                   (2ml)

3M KCl	(666.67ul)
0.25M EDTA	(80ul)
10% NP40	(8ml)
dd H2O	(189.25ml)

**Extraction Buffer (200ml total)**

1M HEPES	(4ml)
3M NaCl	(26.66ml)
0.25M EDTA	(800ul)
dd H2O	(168.54ml)