

Preparation of labeled RNA probes:

I think one key to making a good fluorescent *in situ* is having a smoking hot probe. With enzyme-based (alkaline phosphatase) detection, a mediocre probe can still give you a nice stain because you just let the reaction go longer. But if you don't have very good probes, you might be disappointed with the results of this protocol. What is a very good probe on the alkaline phosphatase reaction-time scale? A probe whose signal is well-developed after 10-20 minutes reaction time, and over-developed by 30 minutes. Truly flaming hot probes are over-developed by 10 minutes. I'm following this fairly standard protocol for making probes.

Preparing and Quantitating DNA Template:

- 1. Linearize 20 μg cDNA plasmid with 5' enzyme (sense strand) in 100 μl . Only 1 μg of linearized template is required in the transcription reaction, so you can digest a lot less. I start with this large digest because I'm making several probes from each template and also like to have a back-up supply.**
- 2. Run out 2 μl on a gel to check for complete cutting**
- 3. Increase volume to 200 μl with TE and add 20 μl 3 M sodium acetate, pH 5.2**
- 4. Perform a phenol/chloroform extraction as you like to do it. Either**
- 5. Extract 1x with phenol and 1x with chloroform OR**
- 6. Extract 2x with phenol/chloroform**

****RNase-free begins here: tips, tubes, reagents, and take a shower before work!**

- 7. Precipitate DNA with 2x volume ethanol**
- 8. Freeze at -20°C . for >30 minutes, then spin at maximum speed for 15 minutes in a microfuge**
- 9. Wash pellet with 200 μl 70% ethanol, spin at maximum speed for 2 minutes, drain completely and air dry**
- 10. Resuspend in 30 μl RNase-free ddH₂O**
- 11. For accurate quantitation by spectrophotometer, dilute 2 μl resuspended template in 100 μl TE, and take the Abs₂₆₀ reading in a microcuvette**

Synthesizing Labeled RNA:

1. Heat template DNA to 55° C. for 2 minutes, then put back on ice
2. Set up this reaction on ice in a RNase-free tube:

| Component | Volume(μ l) |
|-------------------------------|------------------|
| ~1 μ g template DNA | X |
| RNase-free ddH ₂ O | Y |
| 10x transcription buffer | 1.5 |
| 10x hapten-U NTP mix | 1.5 |
| RNase inhibitor | 1.0 |
| RNA polymerase(T7, T3) | 1.5 |
| Total | 15 |

**** Note that the reaction recipe on the Roche product information sheets calls for greater volume of reagents, but the same amount of template: 20 μ l total, instead of 15 μ l as shown here; many protocols call for 10 μ l.**

3. Mix thoroughly with p20 and incubate at 37° C. for 2-2.5 hours.
4. Add 11 μ l RNase-free ddH₂O
5. Take 1 μ l out to put on a gel, then store the reaction at -20° C. or proceed to the fragmentation step
6. Check the reaction product by running it on a regular 0.9% TAE/ethidium bromide gel with a RNA marker. Your RNA product should run in a tight band at about the predicted size and should be at least 10-fold stronger than the DNA template band.

Fragmentation and Precipitation:

1. Add 25 μ l 2x carbonate buffer
2. Mix and incubate at 65° C. for 20-40 minutes (vary to control average probe fragment size)
3. Add 50 μ l stop solution ... optional: take out 3-5 μ l to see the fragmented probe sizes on a gel
4. Add 10 μ l 4 M lithium chloride
5. Add 5 μ l 20 mg/ml tRNA (phenol/chloroform extracted, ethanol precipitated)
6. Add 300 μ l ethanol
7. Vortex and freeze at -20° C. for >30 minutes
8. Spin at maximum speed for 20 minutes in a 4° C. microfuge
9. Wash pellet with 300 μ l 70% ethanol, spin at maximum speed for 2 minutes, drain completely and air dry
10. Dissolve pellet in 200 μ l hybridization solution. Don't let the probe pellet dry for too long, it might become hard to resuspend. Let the probe dissolve on ice for a while, then mix it thoroughly by pipetting with a p200 and vortexing. Probe stocks should be stored at -20° C.